

# K-edge subtraction X-ray imaging with a pixellated spectroscopic detector

Silvia Pani

Department of Physics, University of Surrey



# Summary

- Hyperspectral imaging
- K-edge subtraction X-ray imaging for mammography
  - Methods
  - Results so far
  - Open problems
  - Perspectives

# What is hyperspectral imaging?


- Any imaging technique allowing to retrieve a spectrum of the radiation detected
- Spectral information can be used for
  - Scatter rejection
  - Dual energy techniques
  - Non-conventional techniques (diffraction, fluorescence)
  - Nuclear Medicine
- So far, hyperspectral X-ray imaging methods have been limited due to the limited availability of pixellated spectroscopic detectors
  - Scanning methods are time consuming
- Recent technological developments are making hyperspectral imaging a reality

## In the beginning was HPGe...

- HPGe has excellent energy resolution ( $<0.5$  keV @ 60 keV)
- BUT:
  - Very expensive
  - Needs cooling  $\rightarrow$  bulky and not suitable for in-field applications
  - Very few examples of position-sensitive devices
- HEXITEC collaboration (EPSRC funded):
  - Collaboration between STFC (Rutherford Appleton Labs), Universities of Manchester, Surrey, Durham, Royal Surrey County Hospital
  - Development of pixellated room temperature spectroscopic detectors (CdTe, CZT)
  - Now in its translational phase


# HEXITEC characteristics

- Currently using room temperature/sub-room temperature CdTe sensors
- Current sensors 80x80 pixels, 250  $\mu\text{m}$  pitch, 1 mm thick
  - Tiled arrays of 2 x 2 sensors are under development
- Operated at -500V



# K-edge subtraction imaging

- Typically used for angiography/angiogenesis studies
- Two images are acquired with energies above and below the K-edge of a contrast agent
- Combination and subtraction of the two removes the background
- Problems with conventional KES
  - Increased dose (two images!)
  - Image registration (the patient moves between the two images)
- This is removed with a spectroscopic detector:
  - The images above and below the K-edge are obtained by integrating the spectrum in appropriate ranges



# Iodine KES with HEXITEC

- Aim: study of angiogenesis around breast tumours
- Different concentrations of Niopam® are used
  - Iodine-based contrast agent
- Iodine K-edge: 33.2 keV
  - Higher energy spectra than in conventional mammography
  - W anode, 3 mm Al filtration, 45-50 kVp
- Plus single-photon counting detector
  - Intrinsically low noise
- Significant dose reduction compared to conventional mammography!

# The simplest approach: logarithmic subtraction

- The assumption is that the attenuation coefficient of the background materials does not vary strongly on the two sides of the K-edge
- by subtracting the logarithm of the image below the K-edge from the logarithm of the image above the K-edge the background is removed

$$\begin{aligned} S(i,j) &= \ln(I_{high}^0 / I_{high}(i,j)) - \ln(I_{low}^0 / I_{low}(i,j)) \\ &= \mu_{high}^I x_I + \mu_{high}^{bg} x_{bg} - (\mu_{low}^I x_I + \mu_{low}^{bg} x_{bg}) \\ &\approx (\mu_{high}^I - \mu_{low}^I) x_I \end{aligned}$$


- Problems:
  - Only works if  $\Delta\mu^{bg}$  is negligible across the K-edge
  - Maximum if  $\Delta\mu^I$  is maximum
  - Integration bands must be close to each other
  - Poor background removal if wide bands are chosen



# Finding the optimum integration ranges

- 3 mm tube filled with undiluted Niopam 150 (150  $\mu\text{g/ml}$  I)
- The optimum integration ranges for the images above and below the K-edge result from a trade-off
  - Wide range  $\rightarrow$  high statistics, BUT low contrast + high structural noise
  - Narrow range  $\rightarrow$  high contrast, BUT poor statistics/high noise
- Investigated via contrast-to-noise ratio;  $\text{CNR} = (I - I_b)/\sigma$


The drop in transmission between the two sides of the K-edge is not fully exploited



Width (keV)	CNR
0.5	48
1	76
2	55
5	32
10	14

# Complex test objects

- “non uniform background” phantoms to verify the effectiveness of subtraction algorithms
- Images taken with a W-anode tube, 2 mm filtration, 45 kVp, 7 $\mu$ A
  - Entrance dose ~270  $\mu$ Gy




“single layer phantom”  
1 cm thick



“multi-layer phantom”  
4 cm thick

# Initial K-edge studies

- “non uniform background” phantoms to verify the effectiveness of subtraction algorithms
- Images taken with a W-anode tube, 2 mm filtration, 45 kVp, 7 $\mu$ A
  - Entrance dose  $\sim$ 270  $\mu$ Gy



“single layer phantom”  
1 cm thick




“multi-layer phantom”  
4 cm thick


# Results – logarithmic subtraction

“one layer” phantom

Above K-edge



Below K-edge



Subtracted




“multilayer” phantom




Limited background removal:  
Outside linear regime

# Increasing contrast and improving background subtraction

- When different pixel gain is not taken into account, conservative integration ranges need to be chosen



- The problem is removed by interpolating the spectra for all calibrated pixels so that peaks/K-edges coincide



## Results – optimised calibration

- The procedure allows the choice of integration bands closer to the physical K-edge of contrast agent
  - Higher intrinsic contrast
  - Lower noise due to improved background removal

non


corrected



corrected



Width (keV)	CNR – uncorr	CNR – corr
1.2	31	54
2	26	41
3	31	36



## “Dual energy algorithm” (Lehmann 1981)

- Any image pixel is seen as a vector
- The vector is projected onto a basis {iodine component, water component} and the two projections are retrieved separately, giving rise to a “iodine equivalent” image and a “water equivalent” image


$$\begin{cases} \ln(I_0/I)_{low} = \left(\frac{\mu}{\rho}\right)_{low}^I (\rho x)_I + \left(\frac{\mu}{\rho}\right)_{low}^{water} (\rho x)_{water} \\ \ln(I_0/I)_{high} = \left(\frac{\mu}{\rho}\right)_{high}^I (\rho x)_I + \left(\frac{\mu}{\rho}\right)_{high}^{water} (\rho x)_{water} \end{cases}$$

- Allows quantitative information ( $\mu x$ )
- No pre-assumptions on the position and width of integration bands
  - Proved much more effective than log-subtraction on conventional KES


# Results – Dual energy algorithm (2 keV)

- Background removal is improved!
- Increased CNR (more uniform background)


Log subtracted



Iodine projection




Water projection








# Dual energy optimisation

- “Breast equivalent” test object (Perspex spheres + oil)
- 50 kVp, 3 mm Al (mean energy 33 keV)
- 2  $\mu$ A, 9 min acquisitions  $\rightarrow$  entrance air dose 75  $\mu$ Gy




# Dual energy optimisation - II



- CNR reaches a maximum at wider band than with log subtraction
  - Better background removal → lower structural noise

# Iodine vs water




- Unlike with the log-subtraction algorithm, contrast increases when increasing the bandwidth (but structural noise increases, too)

log subtraction



iodine

water

2 keV




17 keV





# What is going on?

- The two base vectors are no longer orthogonal
- A more accurate treatment will involve changing the system of linear equations for dual energy into a system of integrals




# Open problems - calibration

- The energy calibration of an 80x80 array must be done automatically
- Typically, the user selects a range for peak search using the spectrum averaged on all pixels
- This may cause a problem with pixels with gain significantly off-average → speckles




# Comparison with conventional imaging

- Pairs of spectra were chosen with average energies below and above the K-edge of iodine, respectively
- CNR was calculated to identify the optimum pair of spectra
  - Low: W anode, 45 kVp, 250  $\mu\text{m}$  Tin filter
  - High: W anode, 50 kVp, 12.5 mm Al filter



iodine



water


Tube $\varnothing$ (mm)	CNR – HEXITEC	CNR – conventional
1	25	15
2	55	40
3	68	62

Entrance dose  
Conventional 1.2 mGy  
HEXITEC 75  $\mu\text{Gy}$


# A large test object

- Mean glandular dose 0.5 mGy  
(in conventional  
mammography 1.5  
mGy/image)


LOW




IODINE



HIGH




WATER



Tube + water component of solution


# Very low dose – MGD $1.5 \mu\text{Gy}$



Full spectrum




iodine



low



water



high



# Conclusions

- A pixellated spectroscopic CdTe detector has proven effective in one-shot K-edge subtraction imaging around the Iodine K-edge
  - Dose 15 times lower than conventional imaging with comparable image quality
  - The optimum energy band to be integrated results from a trade-off between signal and noise (statistical and structural)
- Open problems:
  - Limited linearity (to be addressed with a systematic study on pulse shaping parameters)
  - Pixel gain spread
- Future work
  - Improving spatial resolution (sub-pixel) with charge-sharing algorithms
  - Design of a “dynamic” test object for uptake-washout measurements
  - New contrast agents
    - Gold nanoparticles?

# Acknowledgments



- SurreyUni/RSCH
  - Filipa Ferreira
  - Nick Henthorn
  - Sarene Saifuddin
  - James Scuffham
  - Paul Sellin
  - Philipp Stratmann

## The HEXITEC team

- Manchester
  - Bob Cernik
  - Simon Jacques
- RAL
  - Paul Seller
  - Matt Veale
  - Matt Wilson